



ORTHO CARE



Arthroscopy

shoulder implants



FROM WHEELS TO HEELS

Keeping You Strong... Keeping You Healthy... Bone Care For The Lifetime...

VISION

Unparalleled growth and an obsession to take the company to the forefront of the Orthopedic Industry in India with a huge presence across the globe.

MISSION

The mission of ORTHO CARE is to create an innovative and sustainable medium for the development of high-impact and affordable orthopedic implants for the masses and developing countries.

VALUES

Quality: *We shall strive to achieve the highest quality standards that are the hallmark of excellence.*

Customer Focus: *All our efforts and resources shall be channelizing towards meeting and exceeding customer expectations.*

Performance Focus: *We shall have consistent achievement of set business and realize opportunities for accelerating growth.*





Care Anchor Screw Single Loaded Without Needle

Sizes and reference codes

Ref. code	Sizes
CASSL-WN-028	Care Anchor Screw Single Loaded Without Needle 2.8mm
CASSL-N-028	Care Anchor Screw Single Loaded With Needle 2.8mm
CASSL-WN-035	Care Anchor Screw Single Loaded Without Needle 3.5mm
CASSL-N-035	Care Anchor Screw Single Loaded With Needle 3.5mm

Ø3.0
mm

DIAMETER

Care Anchor Peek Screw Single Loaded Without Needle

Sizes and reference codes

Sizes

Care Anchor Peek Screw Single Loaded Without Needle 2.8mm

Care Anchor Peek Screw Double Loaded Without Needle 3.0mm

Care Anchor Peek Screw Double Loaded Without Needle 3.5mm

Ref. code

CAPSSL-WN-028

CAPSDL-WN-030

CAPSDL-WN-035

Ø2.8
mm

DIAMETER

Care Anchor Peek Screw Double Loaded MINI Knotless


Sizes and reference codes

Sizes

Care Anchor Peek Screw Double Loaded Knotless with two piece fiber wires 2.8mm

Ref. code

CAPSSL-DN-028



Care Anchor Screw Double Loaded Without Needle

Sizes and reference codes

Sizes

Care Anchor Screw Double Loaded Without Needle 5.0mm

Care Anchor Screw Double Loaded With Needle 5.0mm

Care Anchor Screw Double Loaded Without Needle 5.5mm

Care Anchor Screw Double Loaded With Needle 5.5mm

Ref. code

CASDL-WN-050

CASDL-N-050

CASDL-WN-055

CASDL-N-055



Care Anchor Peek Screw Double Loaded Without Needle

Sizes and reference codes

Sizes

Care Anchor Peek Screw Single Loaded Without Needle 2.8mm
 Care Anchor Peek Screw Double Loaded Without Needle 3.0mm
 Care Anchor Peek Screw Double Loaded Without Needle 3.5mm
 Care Anchor Peek Screw Double Loaded Without Needle 5.0mm
 Care Anchor Peek Screw Double Loaded Without Needle 5.5mm
 Care Anchor Peek Screw Triple Loaded Without Needle 5.5mm

Ref. code

CAPSSL-WN-028
 CAPSDL-WN-030
 CAPSDL-WN-035
 CAPSDL-WN-050
 CAPSDL-WN-055
 CAPSTL-WN-055



Care Anchor Screw PEEK Knotless


Sizes and reference codes

Sizes

Care Anchor Screw PEEK 5.5mm (Titanium Tip)
 Care Anchor Screw PEEK 5.5mm (Peek Tip)

Ref. code

CAPS-TT-055
 CAPS-PT-055




Care Fiber Wire

Sizes and reference codes

Sizes	Ref. code
Care Fiber Wire White 90cm	CFW-090
Care Fiber Wire White 90cm with Single Needle	CFWSN-090
Care Fiber Wire White with Blue 90cm with Single Sided Needle	CFWBS-090N
Care Fiber Wire White with Blue 90cm with Double Sided Needle	CFWBD-090DN
Care Fiber 100cm with Straight Two Sided Needle 255mm	CFW-100255DN
Care Fiber Wire White with Black 90cm with Single Sided Needle	CFWWB-090SN
Care Fiber Wire White with Black 90cm with Double Sided Needle	CFWWB-090DSN

www.orthocare.in
www.orthoimplants.com




FROM **WHEELS** TO **HEELS**

ORTHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi – 110015, INDIA
Phone : +91-11-41424511 | 41424522 | 45092044 Fax : +91-11-47510708
E-mail : orthocareindia2009@gmail.com Web : www.orthocare.in | www.orthoimplants.com
Factory : A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

USFDA 510K CLEARED PRODUCTS
CE 1023 CERTIFIED PRODUCTS
ISO 13485:2012 CERTIFICATION
INDIAN FDA LICENSE NO. 05/SCP/2011



ORTHO CARE

ACL - PCL Instruments set



Arthroscopy

ACL - PCL Instruments



FROM WHEELS TO HEELS

Keeping You Strong... Keeping You Healthy... Bone Care For The Lifetime...

VISION

Unparalleled growth and an obsession to take the company to the forefront of the Orthopedic Industry in India with a huge presence across the globe.

MISSION


The mission of ORTHO CARE is to create an innovative and sustainable medium for the development of high-impact and affordable orthopedic implants for the masses and developing countries.

VALUES

Quality: *We shall strive to achieve the highest quality standards that are the hallmark of excellence.*

Customer Focus: *All our efforts and resources shall be channelizing towards meeting and exceeding customer expectations.*

Performance Focus: *We shall have consistent achievement of set business and realize opportunities for accelerating growth.*



ACL / PCL reconstruction drill guide system

Sizes and reference codes

Name	Ref. code
- Drill guide	ACPC.01
- Elbow aimer	ACPC.02
- Pcl femoral	ACPC.03
- Pcl tibial	ACPC.04
- Angled bullet	ACPC.05



Drill bit graft passing pin / Beath pin

Sizes and reference codes

Name	Ref. code
- Drill bit	ACPC.06
- Beath pin	ACPC.07



Cannulated tibial reamers

Sizes and reference codes

Size	Ref. code
- 6.0mm , 6.5mm	ACPC.08, ACPC.09
- 7.0mm , 7.5mm	ACPC.10, ACPC.11
- 8.0mm , 8.5mm	ACPC.12, ACPC.13
- 9.0mm , 9.5mm	ACPC.14, ACPC.15
- 10.0mm , 10.5mm	ACPC.16, ACPC.17
- 11.0mm , 11.5mm	ACPC.18, ACPC.19
- 12.0mm	ACPC.20



Cannulated femoral flowertip reamers

Sizes and reference codes

Size	Ref. code
- 6.0mm , 6.5mm	ACPC.21, ACPC.22
- 7.0mm , 7.5mm	ACPC.23, ACPC.24
- 8.0mm , 8.5mm	ACPC.25, ACPC.26
- 9.0mm , 9.5mm	ACPC.27, ACPC.28
- 10.0mm , 10.5mm	ACPC.29, ACPC.30
- 11.0mm , 11.5mm	ACPC.31, ACPC.32
- 12.0mm	ACPC.33



Cannulated endoscopic reamer 4.5mm

Sizes and reference codes

Size	Ref. code
- 4.5mm	ACPC.34



Depth gauge for ACL/PCL reconstruction

Sizes and reference codes


Name	Ref. code
- Depth Gauge	ACPC.35



Femoral aimer

Sizes and reference codes

Size	Ref. code
- 5.0mm	ACPC.36
- 7.0mm	ACPC.37



Tunnel plugs

Sizes and reference codes

Size	Ref. code
- 7-8mm	ACPC.38
- 8-9mm	ACPC.39
- 9-10mm	ACPC.40



Tunnel dialators and T-handle for Tunnel dialators

Sizes and reference codes

Size	Ref. code
- 7.0mm	ACPC.41
- 8.0mm	ACPC.42
- 9.0mm	ACPC.43
- T-handle for Tunnel dialators	ACPC.44



Cannulated screw driver for interference screw

Sizes and reference codes

Name	Ref. code
- Screw driver	ACPC.45



Arthroscopic probes

Sizes and reference codes

Name	Ref. code
- Probe	ACPC.46



Microfracture awls

Sizes and reference codes

Size	Ref. code
- 30°	ACPC.47
- 45°	ACPC.48
- 90°	ACPC.49



Graft prep station

Sizes and reference codes

Name	Ref. code
- Graft board	ACPC.50



Graft sizing block

Sizes and reference codes

Size	Ref. code
- 6-7-8-9-10-11-12	ACPC.51
- 6.5-7.5-8.5-9.5-10.5-11.5-12.5	ACPC.52



Tendon strippers, open / closed

Sizes and reference codes

Size	Ref. code
- 6,0mm open	ACPC.53
- 7,0mm open	ACPC.54
- 6,0mm end	ACPC.55
- 7,0mm end	ACPC.56



Straight / Convex rasp

Sizes and reference codes

Size	Ref. code
- straight	ACPC.57
- convex	ACPC.58



Femoral PCL elevator

Sizes and reference codes

Name	Ref. code
- PCL elevator	ACPC.59



Impactor for spiked ligament staple

Sizes and reference codes

Name	Ref. code
- Impactor	ACPC.60






Graphic case for knee arthroscopy Instruments

Sizes and reference codes

Name	Ref. code
- Graphic case	ACPC.61

www.orthocare.in
www.orthoimplants.com



FROM WHEELS TO HEELS

ORTHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi - 110015, INDIA
Phone : +91- 9810038246 | +91-11-41424511 | 41424522
E-mail : orthocareindia2009@gmail.com | **Web :** www.orthocare.in
Factory : A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

USFDA 510K CLEARED PRODUCTS
CE 1023 CERTIFIED PRODUCTS
ISO 13485 : 2016 CERTIFICATION
LICENCE NO. MFG/MD/2021/000347



ORTHO CARE



Arthroscopy

knee implants



FROM WHEELS TO HEELS

Keeping You Strong... Keeping You Healthy... Bone Care For The Lifetime...

VISION

Unparalleled growth and an obsession to take the company to the forefront of the Orthopedic Industry in India with a huge presence across the globe.

MISSION


The mission of ORTHO CARE is to create an innovative and sustainable medium for the development of high-impact and affordable orthopedic implants for the masses and developing countries.

VALUES

Quality: We shall strive to achieve the highest quality standards that are the hallmark of excellence.

Customer Focus: All our efforts and resources shall be channelizing towards meeting and exceeding customer expectations.

Performance Focus: We shall have consistent achievement of set business and realize opportunities for accelerating growth.



Care Inteference Screw



Sizes and reference codes

		Ref. code
5.0 mm	20mm to 35mm	CIFS-0520 to CIFS-0535
6.0mm	20mm to 35mm	CIFS-0620 to CIFS-0635
7.0mm	20mm to 35mm	CIFS-0720 to CIFS-0735
8.0mm	20mm to 35mm	CIFS-0820 to CIFS-0835
9.0mm	20mm to 35mm	CIFS-0920 to CIFS-0935
10.0mm	20mm to 35mm	CIFS-1020 to CIFS-1035
11.0mm	20mm to 35mm	CIFS-1120 to CIFS-1135

Details

- Specifically designed for endoscopic femoral fixation where graft protection is needed.
- Polished, fully radiused threads provide maximum graft protection.
- Wide-bore cannulation for easy guide wire removal.
- The Taper-Tip style offers sharper initiating threads and a more pronounced taper for faster, easier insertion.
- Designed to reduce the possibility of bone plug migration during placement.
- Round head reduces possibility of graft damage during or after implantation.
- Universal screw for either tibial or femoral applications. Standard 4.5mm hex for easy insertion and removal.
- Cannulated for use with a guide wire.
- Material - Titanium alloy

Care Inteference Screw Round Head



Sizes and reference codes

		Ref. code
7.0mm	20mm to 35mm	CIFSR-0720 to CIFSR-0735
8.0mm	20mm to 35mm	CIFSR-0820 to CIFSR-0835
9.0mm	20mm to 35mm	CIFSR-0920 to CIFSR-0935
10.0mm	20mm to 35mm	CIFSR-1020 to CIFSR-1035
11.0mm	20mm to 35mm	CIFSR-1120 to CIFSR-1135



Details

- Specifically designed for endoscopic femoral fixation where graft protection is needed.
- Polished, fully radiused threads provide maximum graft protection.
- Wide-bore cannulation for easy guide wire removal.
- The Taper-Tip style offers sharper initiating threads and a more pronounced taper for faster, easier insertion.
- Designed to reduce the possibility of bone plug migration during placement.
- Round head reduces possibility of graft damage during or after implantation.
- Universal screw for either tibial or femoral applications. Standard 4.5mm hex for easy insertion and removal.
- Cannulated for use with a guide wire.
- Material - Titanium alloy

Care Inteference Peek Screw



Sizes and reference codes

		Ref. code
5.0 mm	20mm to 35mm	CIFSP-0520 to CIFSP-0535
6.0mm	20mm to 35mm	CIFSP-0620 to CIFSP-0635
7.0mm	20mm to 35mm	CIFSP-0720 to CIFSP-0735
8.0mm	20mm to 35mm	CIFSP-0820 to CIFSP-0835
9.0mm	20mm to 35mm	CIFSP-0920 to CIFSP-0935
10.0mm	20mm to 35mm	CIFSP-1020 to CIFSP-1035
11.0mm	20mm to 35mm	CIFSP-1120 to CIFSP-1135

Details

- Interference Screws provide strong mechanical fixation for both bone patellar tendon-bone (BTB) and soft tissue grafts in ACL and PCL reconstruction.
- Interference Screws feature a thread pattern that allows for a simple surgical technique, with minimal tunnel preparation.
- Non-absorbable interference screw incorporates PEEK polymer from Evonik.
- PEEK material exhibits an ideal combination of strength, stiffness and toughness which surpasses absorbable materials.
- PEEK (Poly-Ether-Ether-Ketone) is a thermoplastic material with excellent biocompatibility and biostability characteristics as per ASTM F2026.
- Interference screws offer the revision ability of an absorbable screw.
- PEEK has inherent radiolucent characteristics that have widely gained confidence and acceptance in the medical community.

Low profile cancellous screw



Sizes and reference codes

		Ref. code
4.0mm	25mm	CIFLPS-0425
	30mm	CIFLPS-0430
	35mm	CIFLPS-0435
	40mm	CIFLPS-0440
6.5mm	25mm	CIFLPS-6525
	30mm	CIFLPS-6530
	35mm	CIFLPS-6535
	40mm	CIFLPS-6540

Details

- The low-profile screw and washer combination provides 31% to 50% lower profile than competitive versions.
- Low profile cancellous screw and spiked washer set provide continuous strength for soft tissue fixation to bone.
- Low profile head to reduce soft tissue irritation.

Spiked Washer For Low Profile Cancellous Screw



Sizes and reference codes

Description	Ref. code
For 4.0mm Cancellous Screw	CSW-0405
For 6.5mm Cancellous Screw	CSW-0657


- Can be used as a post screw for suture or soft tissue fixation.
- Ideal for ACL, PCL, collateral ligament or posteriolateral reconstruction.
- Unique washer design with spikes on the inferior surface provides excellent fixation of soft tissue to bone.
- Material - Titanium alloy.

Suture Disc Slotted



	Ref. code
8mm	CSSD-08
10mm	CSSD-10
12mm	CSSD-12
14mm	CSSD-14

Suture Disc Conical Hole



	Ref. code
8mm	CSCD-08
10mm	CSCD-10
12mm	CSCD-12
14mm	CSCD-14

Care Ligamnet U Shape Stapel



Sizes and reference codes



Ref. code

6.0mm	CLUS-06
8.0mm	CLUS-08
10mm	CLUS-10

Details

- The staple fixation system provides a simple means of attaching tissue to bone.
- Low-profile, wide staple base provides better load distribution.
- Spike-post leg design encourages uninterrupted vascular flow to the underlying tissue.
- Sharp leg points for easier penetration into cortical bone.
- Two pronged implant for superior pull out strength.
- Material-Titanium alloy.

Care Suture Post Bolt



Sizes and reference codes



Ref. code

6.5mm	25mm	CSPB-6525
	30mm	CSPB-6530
	35mm	CSPB-6535

Details

- Used to securely trap attached sutures and bulky suture knots atraumatically under the low-profile head to eliminate soft tissue irritation
- Material-Titanium alloy.

Care Transfix Screw



Ref. code : CTS-01

Details

- The TransFix Screw is rotated into position to securely fixate the graft.
- Self-cutting, self-tapping screw threads allow for easy insertion or removal during revision procedures.
- The TransFix Screw is made from titanium for maximum strength.

CARE Endo Button



Ref. code : CEB-000

Details

- The titanium button is the original endoscopic femoral soft tissue fixation device.
- It can be used for single inclusive incision, soft tissue, or bone-tendon-bone fixation.

CARE Endo Button with loop



Details

- Solid fixation for either primary or auxiliary fixation in ACL or PCL reconstruction.
- Continuous loop of suture braid eliminates the need for knot tying.
- Ideal for a post or screw and washer technique.
- Distal fixation allows circumferential healing of the ACL graft to tunnel wall.
- Significantly enhances soft tissue apposition to bone to speed healing and strengthen the repair.

Sizes and reference codes

Loop size	Ref. code
10mm	CEBL-010
12mm	CEBL-012
15mm	CEBL-015
20mm	CEBL-020
25mm	CEBL-025

ACL/ PCL Reconstruction Set


Ref. code : ACPC-1300.00

Set Items

- 1 ACL / PCL reconstruction drill guide system
- 2 Graft passing 2.4mm Drill bit
- 3 Graft passing pin (Beath pin)
- 4 Cannulated Tibial Reamer,
(Dia. : 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm ,11mm, 12mm)
- 5 Cannulated Femoral Flowertip Reamer
(Dia. : 6.0mm, 7.0mm, 8.0mm, 9.0mm, 10mm ,11mm, 12mm)
- 6 Cannulated Endoscopic Reamer,
(Dia.4.5mm)
- 7 Depth Gauge for ACL/PCL reconstruction
- 8 Femoral Aimer (5mm, 7mm)
- 9 Tunnel Plug (7-8mm, 8-9mm, 9-10mm)
- 10 Tunnel Dilator, (7.0mm, 8.0mm, 9.0mm)
- 11 T-Handle For Tunnel Dilators
- 12 Cannulated Screw Driver For Bio-Lock Interference Screws
- 13 Impactor For Spiked Ligament Staple
- 14 Arthroscopic Probe
- 15 Micro fracture Awl, (30°, 45°, 90°)
- 16 Graft Prep Station
- 17 Graft Sizing Block
- 18 Tendon Stripper, Closed End,6.0mm
- 19 Tendon Stripper, Open End,7.0mm
- 20 Straight Rasp
- 21 Convex Rasp
- 22 Femoral PCL Elevator
- 23 Graphic Case For Knee Arthroscopy Instruments



www.orthocare.in
www.orthoimplants.com



FROM **WHEELS** TO **HEELS**

ORTHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi – 110015, INDIA
Phone : +91-11-41424511 | 41424522 | 45092044 Fax : +91-11-47510708
E-mail : orthocareindia2009@gmail.com Web : www.orthocare.in | www.orthoimplants.com
Factory : A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

USFDA 510K CLEARED PRODUCTS
CE 1023 CERTIFIED PRODUCTS
ISO 13485:2012 CERTIFICATION
INDIAN FDA LICENSE NO. 05/SCP/2011