

CARE DAKSH PFNA NAIL

Implant & Instrument Set



{Deluxe Quality}

CARE DAKSH PFNA NAIL

CARE DAKSH PFNA Nail Implants



DAKSH Proximal Femoral Nail Antirotation (PFNA)
Dia 9mm to 12mm, Length 170mm to 240mm



DAKSH Proximal Femoral Nail Antirotation (PFNA) - Blade
Dia 10.5mm, Length 70mm to 120mm



DAKSH Proximal Femoral Nail Antirotation (PFNA) - Long Right & Left
Dia 9mm to 12mm, Length 320mm to 440mm

(Option - I)



End Cap with Set Screw for PFNA Nail
Length 00mm to 15mm

(Option - II)



End Cap for PFNA Nail
Length 00mm to 15mm



Set Screw for PFNA Nail
Length 00mm



DAKSH PFNA Lag Screw
Length 70mm to 120mm



End Cap for PFNA Lag Screw
Length 00mm to 15mm

Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

CARE DAKSH PFNA Nail Instruments

TRAY-4



CDPFNA-6600.01 **Distal Aiming Device for 170mm & 240mm (Dynamic)**



CDPFNA-6600.02 **Distal Aiming Device for 240mm (Dynamic & Static)**



CDPFNA-6600.03 **Protection Sleeve 16 mm/11 mm & Guide wire sleeve 11 mm /3.2 mm For PFNA Blade**



CDPFNA-6600.04 **Distal Aiming Arm (extension nail)**



CDPFNA-6600.05 **Driving Cap with thread, for Insertion Handle**



CDPFNA-6600.06 **Proximal Aiming Device**



CDPFNA-6600.07 **Distal Locking Guide Bar Part (1)**



CDPFNA-6600.08 **Distal Locking Guide Bar Part (2)**



CDPFNA-6600.09 **PFNA Blade Impactor (Spiral)**



CDPFNA-6600.10 **Hexagonal Screw Driver Tip 3.5mm**



CDPFNA-6600.11 **Insertion Handle**



Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

CARE DAKSH PFNA Nail Instruments

CDPFNA-6600.12

Guide Wire Beaded ϕ 2.5mm




CDPFNA-6600.13

Slide Hammer



TRAY-1



CDPFNA-6600.14

Spanner



CDPFNA-6600.15

Guide wire holder



CDPFNA-6600.16

U block



CDPFNA-6600.17

T-Handle Quick Handle



CDPFNA-6600.18

Trocar 3.2mm



CDPFNA-6600.19

Drill Sleeve ϕ 5.2mm



CDPFNA-6600.20

L-Key SW 3.5mm



CDPFNA-6600.21

L-Key SW 5.0mm



CDPFNA-6600.22

Screw Driver For Nail Connecting Bolt (Small)



Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

CARE DAKSH PFNA Nail Instruments

CDPFNA-6600.23 T-Bit T-Universal wrench

CDPFNA-6600.24 Connecting screws

CDPFNA-6600.25 Connecting screws

CDPFNA-6600.26 Connecting screws

CDPFNA-6600.27 Nail Connecting Bolt

CDPFNA-6600.28 Guide Wire $\phi 3.2\text{mm}$



TRAY-2



CDPFNA-6600.29 Nail extractor

CDPFNA-6600.30 Reamer Cannulated $\phi 11\text{ mm}$ For Trochanteric Cortex

CDPFNA-6600.31 Proximal Reamer Dia 16.5mm

CDPFNA-6600.32 Depth Gauge

CDPFNA-6600.33 Reamer For PFNA Blade 11 mm

CDPFNA-6600.34 Protection Sleeve 8mm, Drill Sleeve 4.0mm, Trocar 4.0mm



CARE DAKSH PFNA Nail Instruments

CDPFNA-6600.35

Spiral Blade Extractor



CDPFNA-6600.36

Needle guide



CDPFNA-6600.37

Measuring scale for spiral blade



CDPFNA-6600.38

Bone Awl



CDPFNA-6600.39

Key for PFNA Blade



CDPFNA-6600.40

Trocar



CDPFNA-6600.41

Protection Sleeve



CDPFNA-6600.42

Rod with thread, for Insertion Handle




CDPFNA-6600.43

Wire Sleeve



TRAY-3




CDPFNA-6600.44

Cleaning Brush



CDPFNA-6600.45

Measuring scale for Nail



CDPFNA-6600.46

Tissue Protector





Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

CARE DAKSH PFNA Nail Instruments

CDPFNA-6600.47	Locator T-Hexagonal Wrench	
CDPFNA-6600.48	Screw Driver For Nail Connecting Bolt	
CDPFNA-6600.49	Drill Bit ø 5.2mm	
CDPFNA-6600.50	Hook Guide	
CDPFNA-6600.51	Clean Needle	
CDPFNA-6600.52	Drill Bit ø 4.0mm	
CDPFNA-6600.53	Cap Hexagonal Wrench, Extension Rod	
CDPFNA-6600.54	Flexible Reamer Dia 9, 10, 11,12mm	

CARE DAKSH PFNA Nail Instrument Set




Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

CARE DAKSH PFNA Nail Instrument Set

Ref. No.	Description	Material	1 SET
CDPFNA – 6600.00	DAKSH PFNA Nail Instrument Set (Deluxe Quality)		
	CONTENTS:-		
	Tray-4		
CDPFNA – 6600.01	Distal Aiming Device for 170mm & 240mm (Dynamic)	AL	1
CDPFNA – 6600.02	Distal Aiming Device for 240mm (Dynamic & Static)	AL	1
CDPFNA – 6600.03	Protection Sleeve 16 mm/11 mm & Guide wire sleeve 11 mm /3.2 mm For PFNA Blade	SS	1 each
CDPFNA – 6600.04	Distal Aiming Arm (extension nail)	SS	1
CDPFNA – 6600.05	Driving Cap with thread, for Insertion Handle	SS	1
CDPFNA – 6600.06	Proximal Aiming Device	AL	1
CDPFNA – 6600.07	Distal Locking Guide Bar Part (1)	SS	1
CDPFNA – 6600.08	Distal Locking Guide Bar Part (2)	SS	1
CDPFNA – 6600.09	PFNA Blade Impactor (Spiral)	SS	1
CDPFNA – 6600.10	Hexagonal Screw Driver Tip 3.5mm	SS	1
CDPFNA – 6600.11	Insertion Handle	SS	1
CDPFNA – 6600.12	Guide Wire Beaded ø2.5mm	SS	1
CDPFNA – 6600.13	Slide Hammer	SS	1
	Tray-1		
CDPFNA – 6600.14	Spanner	SS	1
CDPFNA – 6600.15	Guide wire holder	SS	1
CDPFNA – 6600.16	U block	SS	1
CDPFNA – 6600.17	T-Handle Quick Handle	SS	1
CDPFNA – 6600.18	Trocar 3.2mm	SS	1
CDPFNA – 6600.19	Drill Sleeve ø5.2mm	SS	1
CDPFNA – 6600.20	L-Key SW 3.5mm	SS	1
CDPFNA – 6600.21	L-Key SW 5.0mm	SS	1
CDPFNA – 6600.22	Screw Driver For Nail Connecting Bolt (Small)	SS	1
CDPFNA – 6600.23	T-Bit T-Universal wrench	SS	1
CDPFNA – 6600.24	Connecting screws	SS	2
CDPFNA – 6600.25	Connecting screws	SS	2
CDPFNA – 6600.26	Connecting screws	SS	2
CDPFNA – 6600.27	Nail Connecting Bolt	SS	2
CDPFNA – 6600.28	Guide Wire ø3.2mm	SS	6
	Tray-2		
CDPFNA – 6600.29	Nail extractor	SS	1
CDPFNA – 6600.30	Reamer Cannulated ø 11 mm For Trochanteric Cortex	SS	1
CDPFNA – 6600.31	Proximal Reamer Dia 16.5mm	SS	1
CDPFNA – 6600.32	Depth Gauge	SS	1
CDPFNA – 6600.33	Reamer For PFNA Blade 11 mm	SS	1
CDPFNA – 6600.34	Protection Sleeve 8 ,Drill Sleeve 4.0, Trocar 4.0mm	SS	2 each
CDPFNA – 6600.35	Spiral Blade Extractor	SS	1
CDPFNA – 6600.36	Needle guide	SS	1
CDPFNA – 6600.37	Measuring scale for spiral blade	SS	1
CDPFNA – 6600.38	Bone Awl	SS	1
CDPFNA – 6600.39	Key for PFNA Blade	SS	1
CDPFNA – 6600.40	Trocar	SS	1
CDPFNA – 6600.41	Protection Sleeve	SS	1
CDPFNA – 6600.42	Rod with thread, for Insertion Handle	SS	1
CDPFNA – 6600.43	Wire Sleeve	SS	1 each
	Tray-3		
CDPFNA – 6600.44	Cleaning Brush	SS	1
CDPFNA – 6600.45	Measuring scale for Nail	SS	1
CDPFNA – 6600.46	Tissue Protector	SS	1
CDPFNA – 6600.47	Locator T-Hexagonal Wrench	SS	1
CDPFNA – 6600.48	Screw Driver For Nail Connecting Bolt	SS	1
CDPFNA – 6600.49	Drill Bit ø 5.2mm	SS	1
CDPFNA – 6600.50	Hook Guide	SS	1
CDPFNA – 6600.51	Clean Needle	SS	1
CDPFNA – 6600.52	Drill Bit ø 4.0mm	SS	2
CDPFNA – 6600.53	Cap Hexagonal Wrench, Extension Rod	SS	1 each
CDPFNA – 6600.54	Flexible Reamer Dia 9, 10, 11,12mm	SS	1 each

Specifications, size, shape, color subject to change without notice.
Picture shape if similar to any other manufacturers is just a co-incident.

www.orthocare.in
www.orthoimplants.com



FROM **WHEELS** TO **HEELS**

ORTHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi – 110015, INDIA
Phone : +91- 9810038246 | +91-11-41424511 | 41424522
E-mail : orthocareindia2009@gmail.com | Web : www.orthocare.in
Factory : A-31, Sector A-4, Tronica City Loni, Ghaziabad (U.P.) INDIA

USFDA 510K CLEARED PRODUCTS
CE 1023 CERTIFIED PRODUCTS
ISO 13485:2016 CERTIFICATION
LICENSE NO. MFG/MD/2021/000347